

Product Information	
Catalog No.	FC11326
Product name	Anti-Human CD48 Antibody (156-4H9), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	156-4H9
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Signaling lymphocytic activation molecule 2, TCT.1, CD48 antigen, B-lymphocyte activation marker BLAST-1, BCM1, CD48, BLAST1, Leukocyte antigen MEM-102, SLAMF2, BCM1 surface antigen, SLAM family member 2
Accession	P09326
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.