

Product Information	
Catalog No.	FC11339
Product name	Anti-Human IL8/CXCL8 Antibody (SAA0376), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0376
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	GCP/IL-8 protein II, MDNCF, IL8/NAP1 form III, IL-8(1-77), Neutrophil-activating factor, GCP/IL-8 protein I, GCP/IL-8 protein VI, Interleukin-8, Chemokine (C-X-C motif) ligand 8, Neutrophil-activating protein 1, Monocyte-derived neutrophil chemotactic factor, (Ser-IL-8)72, MDNCF-b, LYAP, GCP/IL-8 protein V, IL8/NAP1 form IV, MONAP, IL8/NAP1 form V, IL8/NAP1 form VI, NAP-1, Granulocyte chemotactic protein 1, GCP/IL-8 protein III, IL8/NAP1 form II, IL8, IL-8, CXCL8, IL8/NAP1 form I, (Ala-IL-8)77, Protein 3-10C, GCP/IL-8 protein IV, Emotakin, GCP-1, C-X-C motif chemokine 8, NAF, T-cell chemotactic factor, MDNCF-c, Monocyte-derived neutrophil-activating peptide, Lymphocyte-derived neutrophil-activating factor
Accession	P10145
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.