

Product Information	
Catalog No.	FC11390
Product name	Anti-Human CD59 Antibody (TS59), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	TS59
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	MSK21, MAC-inhibitory protein, MEM43 antigen, MAC-IP, MIN2, MIRL, HRF20, Protectin, 1F5 antigen, MIC11, MIN3, Membrane inhibitor of reactive lysis, MACIF, 20 kDa homologous restriction factor, CD59, MIN1, CD59 glycoprotein, HRF-20, Membrane attack complex inhibition factor
Accession	P13987
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.