

Product Information	
Catalog No.	FC11450
Product name	Anti-Human TSHR/LGR3 Antibody (3BD10), FITC
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	3BD10
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Thyroid-stimulating hormone receptor, TSH-R, TSHR, Thyrotropin receptor, LGR3
Accession	P16473
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.