

Product Information	
Catalog No.	FC11506
Product name	Anti-Mouse CD73/NT5E Antibody (SAA0267), FITC
Applications	FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0267
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD73, 5'-nucleotidase, NT5E, NT5, 5'-NT, Ecto-5'-nucleotidase, NTE
Accession	Q61503
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.