

Product Information	
Catalog No.	FC11546
Product name	Anti-Human GUCY2C Antibody (SAA0124), FITC
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0124
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	STAR, Intestinal guanylate cyclase, GC-C, GUC2C, STA receptor, Heat-stable enterotoxin receptor, Guanylyl cyclase C, hSTAR, GUCY2C
Accession	P25092
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.