

| Product Information   |                                                                                                                                                                                                                                   |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC11550                                                                                                                                                                                                                           |
| Product name          | Anti-Human CD40/TNFRSF5 Antibody (4D11), FITC                                                                                                                                                                                     |
| Applications          | FCM                                                                                                                                                                                                                               |
| Species reactivity    | Human                                                                                                                                                                                                                             |
| Host species          | Human                                                                                                                                                                                                                             |
| Clonality             | Monoclonal                                                                                                                                                                                                                        |
| Isotype               | IgG4, kappa                                                                                                                                                                                                                       |
| Clone ID              | 4D11                                                                                                                                                                                                                              |
| Form                  | Liquid                                                                                                                                                                                                                            |
| Storage buffer        | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.                                                                                                                                                                                   |
| Concentration         | 0.5 mg/ml                                                                                                                                                                                                                         |
| Target                | Tumor necrosis factor receptor superfamily member 5, B-cell surface antigen CD40, Bp50, CD40L receptor, TNFRSF5, CDw40, CD40                                                                                                      |
| Accession             | P25942                                                                                                                                                                                                                            |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                                                                                                 |
| Stability and Storage | Store at 4°C for 12 months. Protect from light. Do not freeze.                                                                                                                                                                    |
| Note                  | For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only. |