

Product Information	
Catalog No.	FC11566
Product name	Anti-Mouse CD26/DPP4 Antibody (K9/9), FITC
Applications	IF, WB
Species reactivity	Mouse, Rat
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	K9/9
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TP103, Dipeptidyl peptidase IV, ADABP, Dipeptidyl peptidase 4, CD26, DPP4, Dipeptidyl peptidase IV soluble form, Dipeptidyl peptidase IV membrane form, DPP IV, ADCP2, Adenosine deaminase complexing protein 2, ADCP-2, T-cell activation antigen CD26
Accession	P28843
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.