

Product Information	
Catalog No.	FC11605
Product name	Anti-Human CD197/CCR7 Antibody (506E15#), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	506E15#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CDw197, EVI1, CMKBR7, CCR7, C-C CKR-7, CD197, CCR-7, Epstein-Barr virus-induced G-protein coupled receptor 1, MIP-3 beta receptor, EBI1, CC-CKR-7, BLR2, EBV-induced G-protein coupled receptor 1, C-C chemokine receptor type 7
Accession	P32248
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.