

Product Information	
Catalog No.	FC11619
Product name	Anti-Human NOS2/iNOS Antibody (SAA0105), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0105
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Nitric oxide synthase, inducible, NOS2A, Inducible NO synthase, HEP-NOS, Hepatocyte NOS, Peptidyl-cysteine S-nitrosylase NOS2, iNOS, NOS2, Inducible NOS, NOS type II
Accession	P35228
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.