

Product Information	
Catalog No.	FC11624
Product name	Anti-Mouse CD309/KDR/VEGFR-2 Antibody (DC101), FITC
Applications	FCM, WB
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	DC101
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	FLK-1, Kinase insert domain receptor, VEGFR2, Fetal liver kinase 1, CD309, Vascular endothelial growth factor receptor 2, Protein-tyrosine kinase receptor flk-1, FLK1, VEGFR-2, KDR
Accession	P35918
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.