

Product Information	
Catalog No.	FC11718
Product name	Anti-Human CD156b/ADAM17 Antibody (1022C3), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	1022C3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TACE, Snake venom-like protease, ADAM17, CD156b, CSVp, ADAM 17, Disintegrin and metalloproteinase domain-containing protein 17, TNF-alpha convertase, TNF-alpha-converting enzyme
Accession	P78536
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.