

Product Information	
Catalog No.	FC11726
Product name	Anti-Human CD125/IL5RA Antibody (SAA0047), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0047
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL-5R subunit alpha, IL-5 receptor subunit alpha, IL5RA, IL-5RA, IL-5R-alpha, Interleukin-5 receptor subunit alpha, CDw125, CD125, IL5R
Accession	Q01344
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.