

Product Information	
Catalog No.	FC11759
Product name	Anti-Human CD167a/DDR1 Antibody (SAA0826), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0826
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Protein-tyrosine kinase RTK-6, Protein-tyrosine kinase 3A, NTRK4, EDDR1, CAK, Tyrosine kinase DDR, Mammary carcinoma kinase 10, TRKE, Epithelial discoidin domain-containing receptor 1, CD167a, TRK E, Epithelial discoidin domain receptor 1, HGK2, CD167 antigen-like family member A, NEP, MCK-10, Tyrosine-protein kinase CAK, PTK3A, DDR1, Discoidin receptor tyrosine kinase, RTK6, Cell adhesion kinase
Accession	Q08345
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.