

Product Information	
Catalog No.	FC11789
Product name	Anti-Mouse F4/80 Antibody (SAA2157), FITC
Applications	ELISA, FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	SAA2157
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Adhesion G protein-coupled receptor E1, Cell surface glycoprotein F4/80, EGF-like module receptor 1, EGF-like module-containing mucin-like hormone receptor-like 1, EMR1 hormone receptor, Adgre1, Emr1, Gpf480
Accession	Q61549
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.