

Product Information	
Catalog No.	FC11832
Product name	Anti-Human CSPG4/HMW-MAA/MCSP/NG2 Antibody (SAA0155), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0155
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Chondroitin sulfate proteoglycan 4, CSPG4, Chondroitin sulfate proteoglycan NG2, MCSP, Melanoma chondroitin sulfate proteoglycan, Melanoma-associated chondroitin sulfate proteoglycan
Accession	Q6UVK1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.