

Product Information	
Catalog No.	FC11861
Product name	Anti-Human CD303/CLEC4C Antibody (122A2), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	122A2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CLEC4C, C-type lectin superfamily member 7, CLECSF11, C-type lectin domain family 4 member C, Dendritic lectin, CD303, DLEC, HECL, Blood dendritic cell antigen 2, BDCA-2, BDCA2, CLECSF7
Accession	Q8WTT0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.