

Product Information	
Catalog No.	FC11873
Product name	Anti-Human CD270/TNFRSF14 Antibody (36H12), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	36H12
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	HVEA, HVEM, Herpes virus entry mediator A, Tumor necrosis factor receptor superfamily member 14, Tumor necrosis factor receptor-like 2, HveA, TNFRSF14, CD270, TR2, Herpesvirus entry mediator A
Accession	Q92956
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.