

Product Information	
Catalog No.	FC11878
Product name	Anti-Human LINGO1 Antibody (SAA0824), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0824
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	LRRN6A, LERN1, Leucine-rich repeat neuronal protein 1, Leucine-rich repeat and immunoglobulin domain-containing protein 1, LINGO1, Leucine-rich repeat neuronal protein 6A, Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1
Accession	Q96FE5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.