

Product Information	
Catalog No.	FC11894
Product name	Anti-Human LGR4 Antibody (SAA0107), FITC
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0107
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	GPR48, Leucine-rich repeat-containing G-protein coupled receptor 4, LGR4, G-protein coupled receptor 48, RHJ50601
Accession	Q9BXB1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.