

Product Information	
Catalog No.	FC11908
Product name	Anti-Human CD209 Antibody (2B8), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2B8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD209, DC-SIGN1, Dendritic cell-specific ICAM-3-grabbing non-integrin 1, CD209 antigen, CLEC4L, DC-SIGN, C-type lectin domain family 4 member L
Accession	Q9NNX6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.