

Product Information	
Catalog No.	FC11935
Product name	Anti-Mouse CD207 Antibody (2G3k), FITC
Applications	ELISA, FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2G3k
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	C-type lectin domain family 4 member K, Langerin, CD207, Clec4k
Accession	Q8VBX4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.