

Product Information	
Catalog No.	FC11991
Product name	Anti-Human Reticulon-4/RTN4 Antibody (SAA0171), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0171
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	RTN-x, Reticulon-4, Neuroendocrine-specific protein, Neuroendocrine-specific protein C homolog, NOGO, Nogo protein, Neurite outgrowth inhibitor, NSP, Reticulon-5, KIAA0886, RTN4, Foocen
Accession	Q9NQC3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.