

Product Information	
Catalog No.	FC12006
Product name	Anti-Human CD112R/PVRIG Antibody (CPA.7.021), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	CPA.7.021
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	C7orf15, CD112R, CD112 receptor receptor-related immunoglobulin domain-containing protein, PVRIG, Transmembrane protein PVRIG
Accession	Q6DKI7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.