

Product Information	
Catalog No.	FC12054
Product name	Anti-Human LAMP5 Antibody (D1), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	D1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Brain and dendritic cell-associated LAMP, LAMP-5, Brain-associated LAMP-like protein, LAMP5, Lysosome-associated membrane glycoprotein 5, BAD-LAMP, Lysosome-associated membrane protein 5, C20orf103
Accession	Q9UJQ1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.