

Product Information	
Catalog No.	FC12071
Product name	Anti-Human CD175/Tn Antibody (SAA1994), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1994
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Tn antigen saccharide, N-acetyl- α -D-galactosamine, GalNAc α 1-O-Ser/Thr/Tyr, CD175, Tn, Tn-(GalNAc α 1-O-Ser/Thr)
Accession	ChEBI: 40356
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 μ g per million cells in 100 μ L volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.