

Product Information	
Catalog No.	FC12074
Product name	Anti-Human BTN2A1/2 Antibody (SAA2137), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA2137
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Butyrophilin subfamily 2 member A1, BTN2A1, BT2.1, BTF1, Butyrophilin subfamily 2 member A2, BTN2A2, BT2.2, BTF2
Accession	Q7KYR7, Q8WVV5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.