

Product Information	
Catalog No.	FC12090
Product name	Anti-Human CLDN5/Claudin-5 Antibody (M11), APC
Applications	ELISA, FCM
Species reactivity	Human, Monkey
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	M11
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	AWAL, Transmembrane protein deleted in VCFS, TMVCF, TMDVCF, Claudin-5, CLDN5
Accession	O00501
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.