

Product Information	
Catalog No.	FC12135
Product name	Anti-Human LYPD3 Antibody (SAA0109), APC
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	SAA0109
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	C4.4A, Matrigel-induced gene C4 protein, GPI-anchored metastasis-associated protein C4.4A homolog, MIG-C4, LYPD3, Ly6/PLAUR domain-containing protein 3
Accession	O95274
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.