

Product Information	
Catalog No.	FC12172
Product name	Anti-Mouse CD4 Antibody (GK1.5), APC
Applications	FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	GK1.5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	T-cell surface glycoprotein CD4, T-cell differentiation antigen L3T4, T-cell surface antigen T4/Leu-3, CD4
Accession	P06332
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.