

Product Information	
Catalog No.	FC12265
Product name	Anti-Human CD23 Antibody (MHM6), APC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	MHM6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Lymphocyte IgE receptor, C-type lectin domain family 4 member J, Low affinity immunoglobulin epsilon Fc receptor, CD23A, BLAST-2, Immunoglobulin E-binding factor, IGEBF, CD23, Fc-epsilon-RII, CLEC4J, FCE2, FCER2
Accession	P06734
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.