

Product Information	
Catalog No.	FC12297
Product name	Anti-Human CD243/P-Glycoprotein Antibody (C219), APC
Applications	FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	C219
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Phospholipid transporter ABCB1, MDR1, CD243, ATP-binding cassette sub-family B member 1, ABCB1, Multidrug resistance protein 1, P-glycoprotein 1, ATP-dependent translocase ABCB1, PGY1, P-Glycoprotein, P Glycoprotein
Accession	P08183
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.