

Product Information	
Catalog No.	FC12404
Product name	Anti-Human DRD2 Antibody (SAA0772), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	SAA0772
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	D(2) dopamine receptor, Dopamine D2 receptor, DRD2
Accession	P14416
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.