

Product Information	
Catalog No.	FC12407
Product name	Anti-Human MMP9 Antibody (SAA0116), APC
Applications	ELISA, FCM
Species reactivity	Human, Rat
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0116
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Matrix metalloproteinase-9, CLG4B, GELB, 92 kDa type IV collagenase, 92 kDa gelatinase, Gelatinase B, MMP9, MMP-9
Accession	P14780
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.