

Product Information	
Catalog No.	FC12561
Product name	Anti-Mouse CD40/TNFRSF5 Antibody (FGK45), APC
Applications	FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	FGK45
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Tumor necrosis factor receptor superfamily member 5, B-cell surface antigen CD40, Bp50, CD40L receptor, TNFRSF5, CDw40, CD40
Accession	P27512
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.