

Product Information	
Catalog No.	FC12593
Product name	Anti-Human IL12A/IL-12 p35/NKSF1 Antibody (SAA0380), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0380
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL-12A, NK cell stimulatory factor chain 1, IL12A, CLMF p35, Interleukin-12 subunit alpha, IL-12 subunit p35, NKSF1, Cytotoxic lymphocyte maturation factor 35 kDa subunit
Accession	P29459
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.