

Product Information	
Catalog No.	FC12596
Product name	Anti-Human CD154/CD40LG/TNFSF5 Antibody (5C8), APC
Applications	FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	5C8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	sCD40L, Tumor necrosis factor ligand superfamily member 5, T-cell antigen Gp39, CD40L, TNFSF5, CD154, CD40-L, TNF-related activation protein, CD40 ligand, CD40LG, TRAP
Accession	P29965
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.