

Product Information	
Catalog No.	FC12606
Product name	Anti-Human CD32b/FCGR2B Antibody (6G08), APC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	6G08
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IGFR2, FCG2, Fc-gamma RII-b, CDw32, FcRII-b, IgG Fc receptor II-b, CD32, Fc-gamma-RIIb, Low affinity immunoglobulin gamma Fc region receptor II-b, FCGR2B
Accession	P31994
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.