

Product Information	
Catalog No.	FC12635
Product name	Anti-Human ACVRL1 Antibody (SAA0129), APC
Applications	FCM, WB
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2, kappa
Clone ID	SAA0129
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Serine/threonine-protein kinase receptor R3, ACVRL1, Activin receptor-like kinase 1, ALK-1, ACVRLK1, ALK1, SKR3, TGF-B superfamily receptor type I, TSR-I
Accession	P37023
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.