

Product Information	
Catalog No.	FC12653
Product name	Anti-Human CD200 Antibody (B7V3V2), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2, kappa
Clone ID	B7V3V2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	OX-2 membrane glycoprotein, MOX2, MOX1, CD200
Accession	P41217
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.