

Product Information	
Catalog No.	FC12721
Product name	Anti-Human LOX-1 Antibody (15C4), APC
Applications	ELISA, FCM
Species reactivity	Human astrovirus-8
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	15C4
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	C-type lectin domain family 8 member A, LOX-1, Ox-LDL receptor 1, Lectin-like oxidized LDL receptor 1, OLR1, Oxidized low-density lipoprotein receptor 1, Lectin-like oxLDL receptor 1, LOX1, CLEC8A, hLOX-1, Lectin-type oxidized LDL receptor 1
Accession	P78380
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.