

Product Information	
Catalog No.	FC12762
Product name	Anti-Human CDw210b/IL10RB Antibody (SAA2106), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-mFc (IgG2a)
Clone ID	SAA2106
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Interleukin-10 receptor subunit beta, IL-10 receptor subunit beta, IL-10R subunit beta, IL-10RB, Cytokine receptor class-II member 4, Cytokine receptor family 2 member 4, CRF2-4, Interleukin-10 receptor subunit 2, IL-10R subunit 2, IL-10R2, CDw210b, IL10RB, CRFB4, D21S58, D21S66
Accession	Q08334
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.