

Product Information	
Catalog No.	FC12769
Product name	Anti-Mouse CD47/MER6 Antibody (A4), APC
Applications	FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH, mFc
Clone ID	A4
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Leukocyte surface antigen CD47, Integrin-associated protein, IAP, CD47
Accession	Q61735
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.