

Product Information	
Catalog No.	FC12785
Product name	Anti-Human SLC39A6 Antibody (SAA1475), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA1475
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	SLC39A6, LIV1, Zrt- and Irt-like protein 6, Zinc transporter ZIP6, ZIP-6, Estrogen-regulated protein LIV-1, Solute carrier family 39 member 6, ZIP6
Accession	Q13433
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.