

Product Information	
Catalog No.	FC12800
Product name	Anti-Human CD213a2/IL13RA2 Antibody (47#), APC
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	47#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL13R, IL-13RA2, IL-13 receptor subunit alpha-2, IL-13R-alpha-2, Interleukin-13-binding protein, CD213a2, IL-13R subunit alpha-2, IL13RA2, Interleukin-13 receptor subunit alpha-2
Accession	Q14627
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.