

Product Information	
Catalog No.	FC12848
Product name	Anti-Human VTCN1/B7-H4 Antibody (SAA0156), APC
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0156
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	V-set domain-containing T-cell activation inhibitor 1, B7h.5, B7 homolog 4, B7-H4, T-cell costimulatory molecule B7x, B7H4, VTCN1, Protein B7S1, Immune costimulatory protein B7-H4
Accession	Q7Z7D3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.