

Product Information	
Catalog No.	FC12908
Product name	Anti-Human CD360/IL21R Antibody (18A5), APC
Applications	FCM
Species reactivity	Human, Mouse, Rat, Monkey
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	18A5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL21R, NILR, Interleukin-21 receptor, IL-21R, IL-21 receptor, CD360, Novel interleukin receptor
Accession	Q9HBE5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.