

Product Information	
Catalog No.	FC12913
Product name	Anti-Mouse CD248 Antibody (G78), APC
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	G78
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Endosialin, TEM1, CD164L1, Tumor endothelial marker 1, CD248
Accession	Q91V98
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.