

Product Information	
Catalog No.	FC12928
Product name	Anti-Human PDGFC Antibody (SAA0450), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0450
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Spinal cord-derived growth factor, PDGFC receptor-binding form, SCDGF, PDGFC latent form, Fallotein, PDGF-C, VEGF-E, Platelet-derived growth factor C, PDGFC
Accession	Q9NRA1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.