

Product Information	
Catalog No.	FC12937
Product name	Anti-Human IL36RN/IL-1F5 Antibody (SAA0398), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0398
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	FIL1D, Interleukin-1 family member 5, IL-36Ra, IL1HY1, Interleukin-36 receptor antagonist protein, IL1RP3, IL-1RP3, IL-1ra homolog 1, Interleukin-1-like protein 1, IL-1-related protein 3, IL1L1, Interleukin-1 HY1, IL-1HY1, Interleukin-1 delta, FIL1 delta, IL-1 delta, IL-1L1, Interleukin-1 receptor antagonist homolog 1, IL36RN, IL1F5, IL-1F5
Accession	Q9UBH0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.